



asma-1™
Electronic asthma monitor



copd-6™
Electronic COPD screener



lung monitor
Smoking cessation aid



lung age
Lung function home monitors

4000 Respiratory Monitors

Variant	PEF	FEV ₁	FEV ₆	FEV ₁ Ratio	FEF 25-75%	FEV 0.5/0.75	Lung Age	Memory	Personal Best	GOLD	PDF Reports	SpV Reports
asma-1™	●	●	○	○	○	○	○	600	●	○	○	○
asma-1 child version	●	●	○	○	○	●	○	○	●	○	○	○
asma-1 usb, s, bt	●	●	○	○	○	○	○	600	●	○	●	●
copd-6™	○	●	●	●	○	○	●	○	○	●	○	○
copd-6 usb, s, bt	○	●	●	●	○	○	●	○	○	●	●	●
lung monitor	○	●	●	●	○	○	○	200	●	○	○	○
lung monitor usb, s, bt	○	●	●	●	●	○	○	200	●	○	●	●
lung age	○	●	○	○	○	○	●	○	○	○	○	○

● Available ○ Not available

usb - USB bt - Bluetooth™ s - Serial

asma-1™ and copd-6™ are trademarks of Vitalograph Ltd.
The Bluetooth combination mark is a registered trademark and is wholly owned by the Bluetooth SIG.

4000 Series