



# alpha III



## Sample Printout from Thermal Printer

### VITALOGRAPH - ALPHA

DATE: 13/ 9/99  
 NAME: .....  
 REFERENCE NO: 1234  
 AGE: 31  
 HEIGHT: 176 CM.  
 SEX: MALE  
 ETHNIC ORIGIN: CAUCASIAN  
 PREDICTION REFERENCE: ERS 95  
 SOFTWARE REF.NO: 61.201 ISSUE 4  
 CALIBRATION UPDATE DATE: 13/ 9/99  
 CALIBRATION CHECK DATE : 13/ 9/99

PRE : 3 TESTS, BEST WITHIN 2%  
 POST: 3 TESTS, BEST WITHIN 3%

	PRE			POST			CHANGE
	PRED	MERS	%	MERS	%	%	
VC	5.22	4.88	93	5.03	96	3	
IVC	5.22	4.93	94	4.90	94	0	
FVC	4.99	4.96	99	5.02	101	1	
FEV1	4.18	4.24	102	4.09	98	-3	
FEV1 %	82	86	4	81	-1	-5	
PEF	577	582	101	654	113	12	
PEF25-75	4.78	4.68	98	3.80	80	-18	



\*\* INTERPRETATION \*\*  
 NORMAL VENTILATORY FUNCTION  
 NO SIGNIFICANT BRONCHODILATOR  
 RESPONSE.  
 CLINICAL IMPROVEMENT FOLLOWING  
 BRONCHODILATOR MAY OCCUR  
 IN THE ABSENCE OF SPIROMETRIC  
 IMPROVEMENT.  
 CURVE VERIFICATION: SIGNED .....

.... from **Vitalograph**<sup>®</sup>  
 the spirometry people

Vitalograph Ltd. • Maids Moreton • Buckingham • MK18 1SW • England  
 Tel: +44 (0) 1280 827110 • Fax: +44 (0) 1280 823302  
 E-mail: sales@vitalograph.co.uk • Internet: www.vitalograph.co.uk  
 Vitalograph GmbH. • Jacobsenweg 12 • 22525 Hamburg • Germany  
 Tel: (040) 547391-0 • Fax: (040) 547391-40  
 E-mail: info@vitalograph.de • Internet: www.vitalograph.de  
 Vitalograph Inc. • 8347 Quivira Road • Lenexa • Kansas 66215 • USA  
 Tel.: (913) 888 4221 • Fax: (913) 888 4259  
 E-mail: vitcs@vitalograph.com • Internet: www.vitalograph.com  
 Vitalograph (Irl) Ltd. • Gort Road Industrial Estate • Ennis • Co. Clare • Ireland  
 Tel.: (065) 6829611 • Fax: (065) 6829289  
 E-mail: sales@vitalograph.ie • Internet: www.vitalograph.co.uk  
**Vitalograph**<sup>®</sup> and Spirotrac<sup>®</sup> are registered trademarks of Vitalograph Ltd.