



COMPACT

... sample printout
from in-built thermal printer

Vitalograph
... the spirometry people

Vitalograph Ltd. • Maids Moreton • Buckingham • MK18 1SW • England
Tel: +44 (0) 1280 827110 • Fax: +44 (0) 1280 823302 • E-mail: sales@vitalograph.co.uk

Vitalograph GmbH. • Jacobsenweg 12 • 22525 Hamburg • Germany
Tel: (040) 547391-0 • Fax: (040) 547391-40 • E-mail: info@vitalograph.de

Vitalograph Inc. • 13310 West 99th Street • Lenexa • Kansas 66215 • USA
Tel.: (913) 888 4221 • Fax: (913) 888 4259 • E-mail: vitcs@vitalograph.com

Vitalograph (Ireland) Ltd. • Gort Road Business Park • Ennis • Co. Clare • Ireland
Tel.: (065) 6864100 • Fax: (065) 6829289 • E-mail: sales@vitalograph.ie

Internet: www.vitalograph.com

Vitalograph, Spirotrac® and COMPACT are registered trademarks of Vitalograph Ltd.

Vitalograph COMPACT Spirometry Report

ID : S001
 NAME : Adora,
 Helen
 SEX: M HEIGHT: 185cm WEIGHT : 90kg
 BIRTH DATE: 1935/JAN/01 AGE : 70




SESSION ID: 20050311095343
 TEST DATE : 2005/MAR/11 TIME : 09:54
 N° ATTEMPTS: 5
 ATS BEST VALUES AT BTPS

FVC REP. : 0.12L PEF REP.: 1%
 FEV1 REP. : 0.11L

ACCURACY CHECK : 2005/MAR/11
 CALIBRATION UPDATE: 2005/JAN/09
 SERVICE DUE: 2006/JAN/12
 VITALOGRAPH COMPACT NO.: 62431

INDEX		NORM	MEAS	%PRED
		PRED	BEST	
VC	L	4.67	4.50	96
FVC	L	4.91	4.60	94
FEV1	L	3.62	4.03	111
PEF	L/s	9.01	11.07	123
F2575	L/s	2.74	4.82	176
FEF25	L/s	7.60	8.78	116
FEF50	L/s	4.49	5.37	120
FEF75	L/s	1.67	2.39	143
FIVC	L	4.91	4.49	91

* BELOW LOWER LIMIT OF NORMALITY (LLN)



SUGGESTED INTERPRETATION:
Normal ventilatory function.

Signature: _____
 Software Reference: 67700 V1.00
 Printed: 2005/MAR/29, 11:15